

— ABLYBIO, Help Your Research



ITIH4 (YD13068) Rabbit mAb

货号: **AYD15666**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	103kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,A-549
细胞定位	Secreted
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is secreted into the blood, where it is cleaved by plasma kallikrein into two smaller forms. Expression of this gene has been detected only in liver, and it seems to be upregulated during surgical trauma. This gene is part of a cluster of similar genes on chromosome 3. Two transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	3700
基因名	ITIH4
Swiss	Q14624 (https://www.uniprot.org/uniprotkb/Q14624/entry)
别名	ITIH4 (YD13068),ITIH4 (YD13068) Rabbit mAb,ITIH4,Inter-alpha-trypsin inhibitor family heavy chain-related protein,Plasma kallikrein sensitive glycoprotein 120,IHRP,ITIH1,PK120

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)