

— ABLYBIO, Help Your Research



# LRPAP1 (YD12963) Rabbit mAb

货号: **AYD15685**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	41kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,THP-1,HepG2
细胞定位	Rough endoplasmic reticulum lumen, Endoplasmic reticulum-Golgi intermediate compartment lumen, Golgi apparatus, cis-Golgi network, Golgi apparatus lumen, Endosome lumen, Cell surface
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a protein that interacts with the low density lipoprotein (LDL) receptor-related protein and facilitates its proper folding and localization by preventing the binding of ligands. Mutations in this gene have been identified in individuals with myopia 23. Alternative splicing results in multiple transcript variants.
基因ID	4043
基因名	LRPAP1
Swiss	P30533 ( <a href="https://www.uniprot.org/uniprotkb/P30533/entry">https://www.uniprot.org/uniprotkb/P30533/entry</a> )
别名	LRPAP1 (YD12963),LRPAP1 (YD12963) Rabbit mAb,LRPAP1,Low density lipoprotein receptor-related protein-associated protein 1,A2MRAP

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)