

— ABLYBIO, Help Your Research



# SAM68 (YD12946) Rabbit mAb

货号: **AYD15688**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	48kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,HeLa,Jurkat,PC-12,Mouse brain,Rat brain
细胞定位	Nucleus, Cytoplasm, Membrane
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the K homology domain-containing, RNA-binding, signal transduction-associated protein family. The encoded protein appears to have many functions and may be involved in a variety of cellular processes, including alternative splicing, cell cycle regulation, RNA 3'-end formation, tumorigenesis, and regulation of human immunodeficiency virus gene expression. Alternative splicing results in multiple transcript variants.
基因ID	10657
基因名	KHDRBS1
Swiss	Q07666 ( <a href="https://www.uniprot.org/uniprotkb/Q07666/entry">https://www.uniprot.org/uniprotkb/Q07666/entry</a> )
别名	SAM68 (YD12946),SAM68 (YD12946) Rabbit mAb,KHDRBS1,GAP-associated tyrosine phosphoprotein p62,Src-associated in mitosis 68 kDa protein,p21 Ras GTPase-activating protein-associated p62,p68,SAM68

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)