

— ABLYBIO, Help Your Research



Vitamin D Binding protein (YD12868) Rabbit mAb

货号: AYD15703

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF FC
推荐浓度	
理论分子量	53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human serum,Mouse serum,THP-1
细胞定位	Secreted
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene belongs to the albumin gene family. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	2638
基因名	GC
Swiss	P02774 (https://www.uniprot.org/uniprotkb/P02774/entry)
别名	Vitamin D Binding protein (YD12868),Vitamin D Binding protein (YD12868) Rabbit mAb,GC,Gc protein-derived macrophage activating factor,Gc-globulin,Group-specific component,Vitamin D-binding protein-macrophage activating factor

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)