

Plastin L (YD12751) Rabbit mAb

货号: **AYD15726**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	70kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,K-562,DU145,Jurkat,Mouse thymus,Mouse spleen
细胞定位	Cytoplasm, cytoskeleton, Cell junction, Cell projection, ruffle membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	<p>Plastins are a family of actin-binding proteins that are conserved throughout eukaryote evolution and expressed in most tissues of higher eukaryotes. In humans, two ubiquitous plastin isoforms (L and T) have been identified. Plastin 1 (otherwise known as Fimbrin) is a third distinct plastin isoform which is specifically expressed at high levels in the small intestine. The L isoform is expressed only in hemopoietic cell lineages, while the T isoform has been found in all other normal cells of solid tissues that have replicative potential (fibroblasts, endothelial cells, epithelial cells, melanocytes, etc.). However, L-plastin has been found in many types of malignant human cells of non-hemopoietic origin suggesting that its expression is induced accompanying tumorigenesis in solid tissues.</p>
基因ID	3936

基因名	LCP1
Swiss	P13796
别名	Plastin L (YD12751)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn