

— ABLYBIO, Help Your Research



# Casein Kinase 2 beta (YD12740) Rabbit mAb

货号: **AYD15730**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	25kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,SW620,MCF7,Mouse spleen,Mouse kidney,Mouse testis,Rat brain
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes the beta subunit of casein kinase II, a ubiquitous protein kinase which regulates metabolic pathways, signal transduction, transcription, translation, and replication. The enzyme is composed of three subunits, alpha, alpha prime and beta, which form a tetrameric holoenzyme. The alpha and alpha prime subunits are catalytic, while the beta subunit serves regulatory functions. The enzyme localizes to the endoplasmic reticulum and the Golgi apparatus. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	1460
基因名	CSNK2B
Swiss	P67870 ( <a href="https://www.uniprot.org/uniprotkb/P67870/entry">https://www.uniprot.org/uniprotkb/P67870/entry</a> )
别名	Casein Kinase 2 beta (YD12740),Casein Kinase 2 beta (YD12740) Rabbit mAb,CSNK2B,Phosvitin,Protein G 5a,CK2N,G5A

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)