

— ABLYBIO, Help Your Research



# Orai1 (YD12715) Rabbit mAb

货号: **AYD15735**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	33kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,MCF7,Mouse lung,Mouse testis,Mouse spleen,Rat lung,Rat liver
细胞定位	Cell membrane, Basolateral cell membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a membrane calcium channel subunit that is activated by the calcium sensor STIM1 when calcium stores are depleted. This type of channel is the primary way for calcium influx into T-cells. Defects in this gene are a cause of immune dysfunction with T-cell inactivation due to calcium entry defect type 1 (IDTICED1).
------	--

基因ID	84876
基因名	ORAI1
Swiss	Q96D31 ( <a href="https://www.uniprot.org/uniprotkb/Q96D31/entry">https://www.uniprot.org/uniprotkb/Q96D31/entry</a> )
别名	Orai1 (YD12715),Orai1 (YD12715) Rabbit mAb,ORAI1,Protein orai-1,Transmembrane protein 142A,CRACM 1

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)