

— ABLYBIO, Help Your Research



# MTAP (YD12480) Rabbit mAb

货号: **AYD15769**

## 产品信息

|       |   |
|-------|---|
| 反应    | Human,Mouse,Rat   |
| 宿主    | Rabbit  |
| 克隆性   | Monoclonal  |
| 预测反应  |   |
| 应用    | WB IHC-P ICC/IF FC IP   |
| 推荐浓度  |   |
| 理论分子量 | 31kDa   |
| 实测分子量 |   |
| 形式    | Liquid  |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated  |
| 阳性对照  | HeLa,HT-29,A-549,NIH/3T3  |
| 细胞定位  | Cytoplasm, Nucleus  |
| 纯化    | 亲和纯化  |

## 抗原信息



|      |  |
|------|--|
| 抗原信息 |  |
|------|--|

## 靶点信息

|      |  |
|------|--|
| 研究背景 | This gene encodes an enzyme that plays a major role in polyamine metabolism and is important for the salvage of both adenine and methionine. The encoded enzyme is deficient in many cancers because this gene and the tumor suppressor p16 gene are co-deleted. Multiple alternatively spliced transcript variants have been described for this gene, but their full-length natures remain unknown. |
|------|--|

|       |  |
|-------|--|
| 基因ID  | 4507   |
| 基因名   | MTAP   |
| Swiss | Q13126 ( <a href="https://www.uniprot.org/uniprotkb/Q13126/entry">https://www.uniprot.org/uniprotkb/Q13126/entry</a> ) |
| 别名    | MTAP (YD12480),MTAP (YD12480) Rabbit mAb,MTAP,5'-methylthioadenosine phosphorylase,MSAP                                |

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)