

— ABLYBIO, Help Your Research



# Cytokeratin 15 (YD12291) Rabbit mAb

货号: **AYD15789**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	49kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	TE-1,HeLa,SW480,OVCAR-3
细胞定位	
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains and are clustered in a region on chromosome 17q21.2.
------	---

基因ID	3866
基因名	KRT15
Swiss	P19012 ( <a href="https://www.uniprot.org/uniprotkb/P19012/entry">https://www.uniprot.org/uniprotkb/P19012/entry</a> )
别名	Cytokeratin 15 (YD12291),Cytokeratin 15 (YD12291) Rabbit mAb,KRT15,Cytokeratin-15,Keratin-15,KRTB

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)