

— ABLYBIO, Help Your Research



# Retinoic Acid Receptor alpha (YD12279) Rabbit mAb

货号: **AYD15793**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,K-562,Raji,THP-1,SH-SY5Y,Mouse testis,Mouse brain,Rat brain
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene represents a nuclear retinoic acid receptor. The encoded protein, retinoic acid receptor alpha, regulates transcription in a ligand-dependent manner. This gene has been implicated in regulation of development, differentiation, apoptosis, granulopoiesis, and transcription of clock genes. Translocations between this locus and several other loci have been associated with acute promyelocytic leukemia. Alternatively spliced transcript variants have been found for this locus.
基因ID	5914
基因名	RARA
Swiss	P10276 ( <a href="https://www.uniprot.org/uniprotkb/P10276/entry">https://www.uniprot.org/uniprotkb/P10276/entry</a> )
别名	Retinoic Acid Receptor alpha (YD12279),Retinoic Acid Receptor alpha (YD12279) Rabbit mAb,RARA,Nuclear receptor subfamily 1 group B member 1,NR1B1

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)