

— ABLYBIO, Help Your Research



Integrin alpha E (YD12171) Rabbit mAb

货号: **AYD15808**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	130kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse spleen,Mouse thymus,Rat testis,Rat thymus
细胞定位	Membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This gene encodes an I-domain-containing alpha integrin that undergoes post-translational cleavage in the extracellular domain, yielding disulfide-linked heavy and light chains. In combination with the beta 7 integrin, this protein forms the E-cadherin binding integrin known as the human mucosal lymphocyte-1 antigen. This protein is preferentially expressed in human intestinal intraepithelial lymphocytes (IEL), and in addition to a role in adhesion, it may serve as an accessory molecule for IEL activation.
基因ID	3682
基因名	ITGAE
Swiss	P38570 (https://www.uniprot.org/uniprotkb/P38570/entry)
别名	Integrin alpha E (YD12171), Integrin alpha E (YD12171) Rabbit mAb, ITGAE, HML-1 antigen, Integrin alpha-IEL, Mucosal lymphocyte 1 antigen

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)