

— ABLYBIO, Help Your Research



c-Kit (YD12167) Rabbit mAb

货号: **AYD15810**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	110kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Cell membrane, Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes the human homolog of the proto-oncogene c-kit. C-kit was first identified as the cellular homolog of the feline sarcoma viral oncogene v-kit. This protein is a type 3 transmembrane receptor for MGF (mast cell growth factor, also known as stem cell factor). Mutations in this gene are associated with gastrointestinal stromal tumors, mast cell disease, acute myelogenous leukemia, and piebaldism. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3815
基因名	KIT
Swiss	P10721 (https://www.uniprot.org/uniprotkb/P10721/entry)
别名	c-Kit (YD12167),c-Kit (YD12167) Rabbit mAb,KIT,Piebald trait protein,Proto-oncogene c-Kit,Tyrosine-protein kinase Kit,p145 c-kit,v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog,SCFR

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)