

— ABLYBIO, Help Your Research



Prostatic Acid Phosphatase (YD12034) Rabbit mAb

货号: **AYD15836**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse thymus,Rat liver,Rat thymus
细胞定位	Secreted, Cell membrane, Lysosome membrane, Nucleus, Cytoplasm, cytosol
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes an enzyme that catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is synthesized under androgen regulation and is secreted by the epithelial cells of the prostate gland. An alternatively spliced transcript variant encoding a longer isoform has been found for this gene. This isoform contains a transmembrane domain and is localized in the plasma membrane-endosomal-lysosomal pathway.
基因ID	55
基因名	ACP3
Swiss	P15309 (https://www.uniprot.org/uniprotkb/P15309/entry)
别名	Prostatic Acid Phosphatase (YD12034),Prostatic Acid Phosphatase (YD12034) Rabbit mAb,ACP3,5'-nucleotidase,Acid phosphatase 3,Ecto-5'-nucleotidase,Protein tyrosine phosphatase ACP3,Thiamine monophosphatase,ACPP

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)