

— ABLYBIO, Help Your Research



# alpha 1 Antitrypsin (YD12015) Rabbit mAb

货号: **AYD15846**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	47kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,Mouse lung
细胞定位	Secreted, Endoplasmic reticulum, extracellular space, extracellular matrix
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息


研究背景	The protein encoded by this gene is secreted and is a serine protease inhibitor whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. Defects in this gene can cause emphysema or liver disease. Several transcript variants encoding the same protein have been found for this gene.
------	--

基因ID	5265
基因名	SERPINA1
Swiss	P01009 ( <a href="https://www.uniprot.org/uniprotkb/P01009/entry">https://www.uniprot.org/uniprotkb/P01009/entry</a> )
别名	alpha 1 Antitrypsin (YD12015),alpha 1 Antitrypsin (YD12015) Rabbit mAb,SERPINA1,Alpha-1 protease inhibitor,Alpha-1-antiproteinase,Serpin A1,AAT,PI

## 产品验证

Human plasma

kDa  
250 -  
150 -  
100 -  
75 -  
50 - ←  
37 -  
25 -  
20 -  
15 -  
10 -



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)