

— ABLYBIO, Help Your Research



# Cytokeratin 5 (YD11992) Rabbit mAb

货号: **AYD15858**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	62kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375,PC-3
细胞定位	Cytoplasm
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in the basal layer of the epidermis with family member KRT14. Mutations in these genes have been associated with a complex of diseases termed epidermolysis bullosa simplex. The type II cytokeratins are clustered in a region of chromosome 12q12-q13.
基因ID	3852
基因名	KRT5
Swiss	P13647 ( <a href="https://www.uniprot.org/uniprotkb/P13647/entry">https://www.uniprot.org/uniprotkb/P13647/entry</a> )
别名	Cytokeratin 5 (YD11992),Cytokeratin 5 (YD11992) Rabbit mAb,KRT5,58 kDa cytokeratin,Cytokeratin-5,Keratin-5,Type-II keratin Kb5

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)