

— ABLYBIO, Help Your Research



MAGEA3 (YD11982) Rabbit mAb

货号: **AYD15862**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse brain,Rat brain
细胞定位	
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita.
基因ID	4102
基因名	MAGEA3
Swiss	P43357 (https://www.uniprot.org/uniprotkb/P43357/entry)
别名	MAGEA3 (YD11982),MAGEA3 (YD11982) Rabbit mAb,MAGEA3,Antigen MZ2-D,Cancer/testis antigen 1.3,MAGE-3 antigen,MAGE3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)