

MAGEA3 (YD11982) Rabbit mAb

货号: **AYD15862**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P ICC/IF IP |
| 推荐浓度 | |
| 理论分子量 | 35kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | A-549,Mouse brain,Rat brain |
| 细胞定位 | |
| 纯化 | |

抗原信息



| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. |
| 基因ID | 4102 |
| 基因名 | MAGEA3 |
| Swiss | P43357 |

| | |
|----|------------------|
| 别名 | MAGEA3 (YD11982) |
|----|------------------|

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn