

SFRP1 (YD11943) Rabbit mAb

货号: **AYD15876**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,293T,SKOV3,Mouse kidney
细胞定位	Secreted
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. Members of this family act as soluble modulators of Wnt signaling; epigenetic silencing of SFRP genes leads to deregulated activation of the Wnt-pathway which is associated with cancer. This gene may also be involved in determining the polarity of photoreceptor cells in the retina.
基因ID	6422
基因名	SFRP1
Swiss	Q8N474

别名	SFRP1 (YD11943)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn