

— ABLYBIO, Help Your Research



# Glutaminase (YD11822) Rabbit mAb

货号: **AYD15899**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	73kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG, NIH/3T3, Mouse brain, Rat kidney
细胞定位	Mitochondrion, Cytoplasm, cytosol, Mitochondrion matrix
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes the K-type mitochondrial glutaminase. The encoded protein is an phosphate-activated amidohydrolase that catalyzes the hydrolysis of glutamine to glutamate and ammonia. This protein is primarily expressed in the brain and kidney plays an essential role in generating energy for metabolism, synthesizing the brain neurotransmitter glutamate and maintaining acid-base balance in the kidney. Alternative splicing results in multiple transcript variants.
基因ID	2744
基因名	GLS
Swiss	O94925 ( <a href="https://www.uniprot.org/uniprotkb/O94925/entry">https://www.uniprot.org/uniprotkb/O94925/entry</a> )
别名	Glutaminase (YD11822),Glutaminase (YD11822) Rabbit mAb,GLS,K-glutaminase,L-glutamine amidohydrolase,GLS1,KIAA0838

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)