

— ABLYBIO, Help Your Research



SUMO2/3 (YD11792) Rabbit mAb

货号: **AYD15908**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P |
| 推荐浓度 | |
| 理论分子量 | 12kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | MCF7 |
| 细胞定位 | Cytoplasm, Nucleus, PML body |
| 纯化 | 亲和纯化 |


抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | This gene encodes a member of the small ubiquitin-related modifier (SUMO) family of eukaryotic proteins. The encoded protein is covalently conjugated to other proteins via a post-translation modification known as sumoylation. Sumoylation may play a role in a wide variety of cellular processes, including nuclear transport, DNA replication and repair, mitosis, transcriptional regulation, and signal transduction. Alternatively spliced transcript variants encoding distinct proteins have been described. |
| 基因ID | 6612 |
| 基因名 | SUMO3 |
| Swiss | P55854 (https://www.uniprot.org/uniprotkb/P55854/entry) |
| 别名 | SUMO2/3 (YD11792),SUMO2/3 (YD11792) Rabbit mAb,SUMO3,SMT3 homolog 1,SUMO-2,Ubiquitin-like protein SMT3A,SMT3A |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)