

— ABLYBIO, Help Your Research



Trypsin (YD11672) Rabbit mAb

货号: **AYD15933**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	27kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Mouse spleen,Mouse kidney
细胞定位	Secreted, extracellular space
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a trypsinogen, which is a member of the trypsin family of serine proteases. This enzyme is secreted by the pancreas and cleaved to its active form in the small intestine. It is active on peptide linkages involving the carboxyl group of lysine or arginine. Mutations in this gene are associated with hereditary pancreatitis. This gene and several other trypsinogen genes are localized to the T cell receptor beta locus on chromosome 7.
基因ID	5644
基因名	PRSS1
Swiss	P07477 (https://www.uniprot.org/uniprotkb/P07477/entry)
别名	Trypsin (YD11672),Trypsin (YD11672) Rabbit mAb,PRSS1,Anionic trypsin I,Anionic trypsin-I,Beta-trypsin,Cationic trypsinogen,Pretrypsinogen I,Trypsin I,Trypsin-1,TRP1,TRY1,TRYP1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)