

— ABLYBIO, Help Your Research



# HNF1 beta (YD11651) Rabbit mAb

货号: **AYD15937**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	61kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Rat liver
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the homeodomain-containing superfamily of transcription factors. The protein binds to DNA as either a homodimer, or a heterodimer with the related protein hepatocyte nuclear factor 1-alpha. The gene has been shown to function in nephron development, and regulates development of the embryonic pancreas. Mutations in this gene result in renal cysts and diabetes syndrome and noninsulin-dependent diabetes mellitus, and expression of this gene is altered in some types of cancer. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	6928
基因名	HNF1B
Swiss	P35680 ( <a href="https://www.uniprot.org/uniprotkb/P35680/entry">https://www.uniprot.org/uniprotkb/P35680/entry</a> )
别名	HNF1 beta (YD11651),HNF1 beta (YD11651) Rabbit mAb,HNF1B,Homeoprotein LFB3,Transcription factor 2,Variant hepatic nuclear factor 1,TCF2

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)