

— ABLYBIO, Help Your Research



IL-8 (YD11649) Rabbit mAb

货号: **AYD15938**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	11kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1,Recombinant Human IL-8 Protein
细胞定位	Secreted
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the CXC chemokine family. This chemokine is one of the major mediators of the inflammatory response. This chemokine is secreted by several cell types. It functions as a chemoattractant, and is also a potent angiogenic factor. This gene is believed to play a role in the pathogenesis of bronchiolitis, a common respiratory tract disease caused by viral infection. This gene and other ten members of the CXC chemokine gene family form a chemokine gene cluster in a region mapped to chromosome 4q.
基因ID	3576
基因名	CXCL8
Swiss	P10145 (https://www.uniprot.org/uniprotkb/P10145/entry)
别名	IL-8 (YD11649),IL-8 (YD11649) Rabbit mAb,CXCL8,C-X-C motif chemokine 8,Chemokine (C-X-C motif) ligand 8,Emoxtakin,Granulocyte chemotactic protein 1,Monocyte-derived neutrophil chemotactic factor,Monocyte-derived neutrophil-activating peptide

产品验证

Recombinant Human IL-8 protein



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)