

— ABLYBIO, Help Your Research



S100A4 (YD11613) Rabbit mAb

货号: **AYD15950**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa cells, A549 cells, A375 cells, PC-3 cells, NIH/3T3 cells
细胞定位	Secreted, Nucleus, Cytoplasm
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in motility, invasion, and tubulin polymerization. Chromosomal rearrangements and altered expression of this gene have been implicated in tumor metastasis. Multiple alternatively spliced variants, encoding the same protein, have been identified.
基因ID	6275
基因名	S100A4
Swiss	P26447 (https://www.uniprot.org/uniprotkb/P26447/entry)
别名	S100A4 (YD11613),S100A4 (YD11613) Rabbit mAb,S100A4,Calvasculin,Metastasin,Placental calcium-binding protein,Protein Mts1,S100 calcium-binding protein A4,CAPL,MTS1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)