

— ABLYBIO, Help Your Research



Fibrinogen alpha chain (YD11541) Rabbit mAb

货号: **AYD15971**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	95kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Rat liver
细胞定位	Secreted
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes the alpha subunit of the coagulation factor fibrinogen, which is a component of the blood clot. Following vascular injury, the encoded preproprotein is proteolytically processed by thrombin during the conversion of fibrinogen to fibrin. Mutations in this gene lead to several disorders, including dysfibrinogenemia, hypofibrinogenemia, afibrinogenemia and renal amyloidosis. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that undergoes proteolytic processing.
基因ID	2243
基因名	FGA
Swiss	P02671 (https://www.uniprot.org/uniprotkb/P02671/entry)
别名	Fibrinogen alpha chain (YD11541),Fibrinogen alpha chain (YD11541) Rabbit mAb,FGA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)