

— ABLYBIO, Help Your Research



Progesterone Receptor (YD11399) Rabbit mAb

货号: **AYD16004**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	99kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Mouse brain,Rat brain
细胞定位	Nucleus, Cytoplasm, Mitochondrion outer membrane
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the steroid receptor superfamily. The encoded protein mediates the physiological effects of progesterone, which plays a central role in reproductive events associated with the establishment and maintenance of pregnancy. This gene uses two distinct promoters and translation start sites in the first exon to produce several transcript variants, both protein coding and non-protein coding. Two of the isoforms (A and B) are identical except for an additional 165 amino acids found in the N-terminus of isoform B and mediate their own response genes and physiologic effects with little overlap.
基因ID	5241
基因名	PGR
Swiss	P06401 (https://www.uniprot.org/uniprotkb/P06401/entry)
别名	Progesterone Receptor (YD11399),Progesterone Receptor (YD11399) Rabbit mAb,PGR,Nuclear receptor subfamily 3 group C member 3,NR3C3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)