

— ABLYBIO, Help Your Research



# Integrin alpha 2 (YD11384) Rabbit mAb

货号: **AYD16006**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	129kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,HT-1080
细胞定位	Membrane
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes the alpha subunit of a transmembrane receptor for collagens and related proteins. The encoded protein forms a heterodimer with a beta subunit and mediates the adhesion of platelets and other cell types to the extracellular matrix. Loss of the encoded protein is associated with bleeding disorder platelet-type 9. Antibodies against this protein are found in several immune disorders, including neonatal alloimmune thrombocytopenia. This gene is located adjacent to a related alpha subunit gene. Alternative splicing results in multiple transcript variants.
基因ID	3673
基因名	ITGA2
Swiss	P17301 ( <a href="https://www.uniprot.org/uniprotkb/P17301/entry">https://www.uniprot.org/uniprotkb/P17301/entry</a> )
别名	Integrin alpha 2 (YD11384),Integrin alpha 2 (YD11384) Rabbit mAb,ITGA2,CD49 antigen-like family member B,Collagen receptor,Platelet membrane glycoprotein Ia,VLA-2 subunit alpha,CD49B

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)