

MALT1 (YD11376) Rabbit mAb

货号: **AYD16007**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	92kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T
细胞定位	Cytoplasm, perinuclear region, Nucleus
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene has been found to be recurrently rearranged in chromosomal translocation with two other genes - baculoviral IAP repeat-containing protein 3 (also known as apoptosis inhibitor 2) and immunoglobulin heavy chain locus - in mucosa-associated lymphoid tissue lymphomas. The protein encoded by this gene may play a role in NF-kappaB activation. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene.
基因ID	10892
基因名	MALT1
Swiss	Q9UDY8

别名	MALT1 (YD11376)
----	-----------------

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn