

— ABLYBIO, Help Your Research



Stathmin 1 (YD11296) Rabbit mAb

货号: **AYD16027**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,C6
细胞定位	Cytoplasm, cytoskeleton
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene belongs to the stathmin family of genes. It encodes a ubiquitous cytosolic phosphoprotein proposed to function as an intracellular relay integrating regulatory signals of the cellular environment. The encoded protein is involved in the regulation of the microtubule filament system by destabilizing microtubules. It prevents assembly and promotes disassembly of microtubules. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3925
基因名	STMN1
Swiss	P16949 (https://www.uniprot.org/uniprotkb/P16949/entry)
别名	Stathmin 1 (YD11296),Stathmin 1 (YD11296) Rabbit mAb,STMN1,Leukemia-associated phosphoprotein p18,Metablastin,Oncoprotein 18,Phosphoprotein p19,Prosolin,Protein Pr22,pp17,C1orf215,LAP18,OP18

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)