

CD48 (YD11263) Rabbit mAb

货号: **AYD16035**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P IP |
| 推荐浓度 | |
| 理论分子量 | 28kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse pancreas,Mouse spleen,Rat pancreas |
| 细胞定位 | Cell membrane, Membrane raft, Secreted |
| 纯化 | |

抗原信息


| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|------|--|
| 研究背景 | This gene encodes a member of the CD2 subfamily of immunoglobulin-like receptors which includes SLAM (signaling lymphocyte activation molecules) proteins. The encoded protein is found on the surface of lymphocytes and other immune cells, dendritic cells and endothelial cells, and participates in activation and differentiation pathways in these cells. The encoded protein does not have a transmembrane domain, however, but is held at the cell surface by a GPI anchor via a C-terminal domain which maybe cleaved to yield a soluble form of the receptor. Multiple transcript variants encoding different isoforms have been found for this gene. |
| 基因ID | 962 |
| 基因名 | CD48 |

| | |
|-------|----------------|
| Swiss | P09326 |
| 别名 | CD48 (YD11263) |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn