

— ABLYBIO, Help Your Research



CD48 (YD11263) Rabbit mAb

货号: **AYD16035**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	28kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse pancreas,Mouse spleen,Rat pancreas
细胞定位	Cell membrane, Membrane raft, Secreted
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the CD2 subfamily of immunoglobulin-like receptors which includes SLAM (signaling lymphocyte activation molecules) proteins. The encoded protein is found on the surface of lymphocytes and other immune cells, dendritic cells and endothelial cells, and participates in activation and differentiation pathways in these cells. The encoded protein does not have a transmembrane domain, however, but is held at the cell surface by a GPI anchor via a C-terminal domain which maybe cleaved to yield a soluble form of the receptor. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	962
基因名	CD48
Swiss	P09326 (https://www.uniprot.org/uniprotkb/P09326/entry)
别名	CD48 (YD11263),CD48 (YD11263) Rabbit mAb,CD48,B-lymphocyte activation marker BLAST-1,BCM1 surface antigen,Leukocyte antigen MEM-102,SLAM family member 2,Signaling lymphocytic activation molecule 2,TCT.1,BCM1,BLAST1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)