

— ABLYBIO, Help Your Research



Hsp70 (YD11253) Rabbit mAb

货号: **AYD16041**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	70kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,BxPC-3,293T,A-549,Mouse brain,Mouse kidney,Mouse lung,Rat lung
细胞定位	Cytoplasm, Nucleus, cytoskeleton, microtubule organizing center, centrosome, Secreted
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This intronless gene encodes a 70kDa heat shock protein which is a member of the heat shock protein 70 family. In conjunction with other heat shock proteins, this protein stabilizes existing proteins against aggregation and mediates the folding of newly translated proteins in the cytosol and in organelles. It is also involved in the ubiquitin-proteasome pathway through interaction with the AU-rich element RNA-binding protein 1. The gene is located in the major histocompatibility complex class III region, in a cluster with two closely related genes which encode similar proteins.
基因ID	3303
基因名	HSPA1A
Swiss	P0DMV8 (https://www.uniprot.org/uniprotkb/P0DMV8/entry)
别名	Hsp70 (YD11253),Hsp70 (YD11253) Rabbit mAb,HSPA1A,Heat shock 70 kDa protein 1,Heat shock protein family A member 1A,HSP72,HSPA1,HSX70

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)