

— ABLYBIO, Help Your Research



HNF-4-alpha (YD11175) Rabbit mAb

货号: **AYD16063**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480
细胞定位	Nucleus
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms.
基因ID	3172
基因名	HNF4A
Swiss	P41235 (https://www.uniprot.org/uniprotkb/P41235/entry)
别名	HNF-4-alpha (YD11175),HNF-4-alpha (YD11175) Rabbit mAb,HNF4A,Nuclear receptor subfamily 2 group A member 1,Transcription factor 14,Transcription factor HNF-4,HNF4,NR2A1,TCF14

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)