

— ABLYBIO, Help Your Research



5T4 (YD11154) Rabbit mAb

货号: **AYD16068**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	46kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse heart,Rat brain,Rat heart,Rat liver
细胞定位	Cell membrane
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a leucine-rich transmembrane glycoprotein that may be involved in cell adhesion. The encoded protein is an oncofetal antigen that is specific to trophoblast cells. In adults this protein is highly expressed in many tumor cells and is associated with poor clinical outcome in numerous cancers. Alternative splicing in the 5' UTR results in multiple transcript variants that encode the same protein.
------	---

基因ID	7162
基因名	TPBG
Swiss	Q13641 (https://www.uniprot.org/uniprotkb/Q13641/entry)
别名	5T4 (YD11154),5T4 (YD11154) Rabbit mAb,TPBG,5T4 oncofetal antigen,5T4 oncofetal trophoblast glycoprotein,M6P1,Wnt-activated inhibitory factor 1,5T4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)