

— ABLYBIO, Help Your Research



# CD8 alpha (YD19949) Rabbit mAb

货号: **AYD16107**

## 产品信息

反应	Human, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC ICC/IF FC
推荐浓度	
理论分子量	26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	THP-1, Mouse spleen
细胞定位	Cell membrane, Secreted
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The CD8 antigen is a cell surface glycoprotein found on most cytotoxic T lymphocytes that mediates efficient cell-cell interactions within the immune system. The CD8 antigen acts as a coreceptor with the T-cell receptor on the T lymphocyte to recognize antigens displayed by an antigen presenting cell in the context of class I MHC molecules. The coreceptor functions as either a homodimer composed of two alpha chains or as a heterodimer composed of one alpha and one beta chain. Both alpha and beta chains share significant homology to immunoglobulin variable light chains. This gene encodes the CD8 alpha chain. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	925
基因名	CD8A
Swiss	P01732 ( <a href="https://www.uniprot.org/uniprotkb/P01732/entry">https://www.uniprot.org/uniprotkb/P01732/entry</a> )
别名	CD8 alpha (YD19949),CD8 alpha (YD19949) Rabbit mAb,CD8A,T-lymphocyte differentiation antigen T8/Leu-2,MAL

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)