

SLC27A4/FATP4 (YD19310) Rabbit mAb

货号: **AYD16126**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	72kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Endoplasmic reticulum membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of a family of fatty acid transport proteins, which are involved in translocation of long-chain fatty acids across the plasma membrane. This protein is expressed at high levels on the apical side of mature enterocytes in the small intestine, and appears to be the principal fatty acid transporter in enterocytes. Clinical studies suggest this gene as a candidate gene for the insulin resistance syndrome. Mutations in this gene have been associated with ichthyosis prematurity syndrome.
基因ID	10999
基因名	SLC27A4
Swiss	Q6P1M0

别名	SLC27A4/FATP4 (YD19310)
----	-------------------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn