

Cyclophilin A (YD35986) Rabbit mAb

货号: **AYD16140**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	18kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,Jurkat,HepG2,HeLa,NIH/3T3,Mouse lung,mouse spleen,rat spleen,mouse brain,rat brain
细胞定位	Cytoplasm, Secreted, Nucleus, Cell membrane
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. The encoded protein is a cyclosporin binding-protein and may play a role in cyclosporin A-mediated immunosuppression. The protein can also interact with several HIV proteins, including p55 gag, Vpr, and capsid protein, and has been shown to be necessary for the formation of infectious HIV virions. Multiple pseudogenes that map to different chromosomes have been reported.
基因ID	5478
基因名	PPIA

Swiss	P62937
别名	Cyclophilin A (YD35986)

产品验证

Mouse spleen

kDa
180 -
150 -
100 -
70 -
55 -
40 -
35 -
25 -



实验步骤

访问官网浏览详情: www.ablybio.cn