

APOF (YD36037) Rabbit mAb

货号: **AYD16272**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Minor apolipoprotein that associates with LDL. Inhibits cholesteryl ester transfer protein (CETP) activity and appears to be an important regulator of cholesterol transport. Also associates to a lesser degree with VLDL, Apo-AI and Apo-AII
基因ID	319
基因名	APOF
Swiss	Q13790
别名	APOF (YD36037)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn