

BRMS1 (YD36041) Rabbit mAb

货号: **AYD16276**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	28kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse testis
细胞定位	Nucleus, Cytoplasm
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene reduces the metastatic potential, but not the tumorigenicity, of human breast cancer and melanoma cell lines. The protein encoded by this gene localizes primarily to the nucleus and is a component of the mSin3a family of histone deacetylase complexes (HDAC). The protein contains two coiled-coil motifs and several imperfect leucine zipper motifs. Alternative splicing results in two transcript variants encoding different isoforms.
基因ID	25855
基因名	BRMS1
Swiss	Q9HCU9

别名	BRMS1 (YD36041)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn