

— ABLYBIO, Help Your Research



CPEB1 (YD36073) Rabbit mAb

货号: **AYD16308**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	63kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,K-562,Mouse brain,Rat liver
细胞定位	Cytoplasm, Nucleus, P-body, Cytoplasmic granule, Synapse, Membrane, Postsynaptic density, Cell project ion, dendrite
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the cytoplasmic polyadenylation element binding protein family. This highly conserved protein binds to a specific RNA sequence, called the cytoplasmic polyadenylation element, found in the 3' untranslated region of some mRNAs. The encoded protein functions in both the cytoplasm and the nucleus. It is involved in the regulation of mRNA translation, as well as processing of the 3' untranslated region, and may play a role in cell proliferation and tumorigenesis. Alternative splicing results in multiple transcript variants.
基因ID	64506
基因名	CPEB1
Swiss	Q9BZB8 (https://www.uniprot.org/uniprotkb/Q9BZB8/entry)
别名	CPEB1 (YD36073),CPEB1 (YD36073) Rabbit mAb,CPEB1,CPEB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)