

CDH6 (YD36147) Rabbit mAb

货号: **AYD16382**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	88kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse brain,Rat kidney,Rat brain
细胞定位	Cell membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis.
基因ID	1004
基因名	CDH6
Swiss	P55285

别名	CDH6 (YD36147)
----	----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn