

— ABLYBIO, Help Your Research



# CDH6 (YD36147) Rabbit mAb

货号: **AYD16382**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	88kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse brain,Rat kidney,Rat brain
细胞定位	Cell membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis.
基因ID	1004
基因名	CDH6
Swiss	P55285 ( <a href="https://www.uniprot.org/uniprotkb/P55285/entry">https://www.uniprot.org/uniprotkb/P55285/entry</a> )
别名	CDH6 (YD36147),CDH6 (YD36147) Rabbit mAb,CDH6,Kidney cadherin

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)