

— ABLYBIO, Help Your Research



# UBE4A (YD36270) Rabbit mAb

货号: **AYD16505**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	123kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the U-box ubiquitin ligase family. The encoded protein is involved in mult iubiquitin chain assembly and plays a critical role in chromosome condensation and separation through t he polyubiquitination of securin. Autoantibodies against the encoded protein may be markers for sclerode rma and Crohn's disease. A pseudogene of this gene is located on the long arm of chromosome 3. Alterna tively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	9354
基因名	UBE4A
Swiss	Q14139 ( <a href="https://www.uniprot.org/uniprotkb/Q14139/entry">https://www.uniprot.org/uniprotkb/Q14139/entry</a> )
别名	UBE4A (YD36270),UBE4A (YD36270) Rabbit mAb,UBE4A,RING-type E3 ubiquitin transferase E4 A,KIAA012 6

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)