

UBE2B (YD36271) Rabbit mAb

货号: **AYD16506**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,SW480,Jurkat,NIH/3T3,Mouse heart,Mouse brain,Mouse skeletal muscle,Rat heart,Rat brain
细胞定位	Cell membrane, Nucleus
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is required for post-replicative DNA damage repair. Its protein sequence is 100% identical to the mouse, rat, and rabbit homologs, which indicates that this enzyme is highly conserved in eukaryotic evolution.
基因ID	7320
基因名	UBE2B

Swiss	P63146
别名	UBE2B (YD36271)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn