

— ABLYBIO, Help Your Research



# VPS28 (YD36276) Rabbit mAb

货号: **AYD16511**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	25kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,U-251MG,293T,HepG2,Mouse brain,Mouse kidney,Mouse liver,Mouse spleen,Rat brain
细胞定位	Cell membrane, Late endosome membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a protein subunit of the ESCRT-I complex (endosomal complexes required for transport), which functions in the transport and sorting of proteins into subcellular vesicles. This complex can also be hijacked to facilitate the budding of enveloped viruses from the cell membrane. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	---

基因ID	51160
基因名	VPS28
Swiss	Q9UK41 ( <a href="https://www.uniprot.org/uniprotkb/Q9UK41/entry">https://www.uniprot.org/uniprotkb/Q9UK41/entry</a> )
别名	VPS28 (YD36276),VPS28 (YD36276) Rabbit mAb,VPS28,ESCRT-I complex subunit VPS28

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)