

— ABLYBIO, Help Your Research



RBM8A (YD36282) Rabbit mAb

货号: **AYD16517**

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,293T,Jurkat,Mouse lung,Mouse brain,Rat lung,Rat brain,Rat spleen
细胞定位	Nucleus, Nucleus speckle, Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein with a conserved RNA-binding motif. The protein is found predominantly in the nucleus, although it is also present in the cytoplasm. It is preferentially associated with mRNAs produced by splicing, including both nuclear mRNAs and newly exported cytoplasmic mRNAs. It is thought that the protein remains associated with spliced mRNAs as a tag to indicate where introns had been present, thus coupling pre- and post-mRNA splicing events. Previously, it was thought that two genes encode this protein, RBM8A and RBM8B; it is now thought that the RBM8B locus is a pseudogene. There are two alternate translation start codons with this gene, which result in two forms of the protein. An allele mutation and a low-frequency noncoding single-nucleotide polymorphism (SNP) in this gene cause thrombocytopenia-absent radius (TAR) syndrome. [provided by RefSeq, Jul 2013]
基因ID	9939
基因名	RBM8A
Swiss	Q9Y5S9 (https://www.uniprot.org/uniprotkb/Q9Y5S9/entry)
别名	RBM8A (YD36282),RBM8A (YD36282) Rabbit mAb,RBM8A,Binder of OVCA1-1,RNA-binding motif protein 8 A,RNA-binding protein Y14,Ribonucleoprotein RBM8A,RBM8

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)