

— ABLYBIO, Help Your Research



COL18A1 (YD35740) Rabbit mAb

货号: **AYD16540**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	178kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted, extracellular space, extracellular matrix, basement membrane
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	Probably plays a major role in determining the retinal structure as well as in the closure of the neural tube
基因ID	80781

基因名	COL18A1
Swiss	P39060 (https://www.uniprot.org/uniprotkb/P39060/entry)
别名	COL18A1 (YD35740),COL18A1 (YD35740) Rabbit mAb,COL18A1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)