

— ABLYBIO, Help Your Research



TIMP3 (YD35940) Rabbit mAb

货号: **AYD16604**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	24kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat kidney
细胞定位	Secreted, extracellular space, extracellular matrix
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene belongs to the TIMP gene family. The proteins encoded by this gene family are inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix (ECM). Expression of this gene is induced in response to mitogenic stimulation and this netrin domain-containing protein is localized to the ECM. Mutations in this gene have been associated with the autosomal dominant disorder Sorsby's fundus dystrophy.
基因ID	7078
基因名	TIMP3
Swiss	P35625 (https://www.uniprot.org/uniprotkb/P35625/entry)
别名	TIMP3 (YD35940),TIMP3 (YD35940) Rabbit mAb,TIMP3,Protein MIG-5,Tissue inhibitor of metalloproteinases 3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)