

TIMP3 (YD35940) Rabbit mAb

货号: **AYD16604**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	24kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat kidney
细胞定位	Secreted, extracellular space, extracellular matrix
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene belongs to the TIMP gene family. The proteins encoded by this gene family are inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix (ECM). Expression of this gene is induced in response to mitogenic stimulation and this netrin domain-containing protein is localized to the ECM. Mutations in this gene have been associated with the autosomal dominant disorder Sorsby's fundus dystrophy.
基因ID	7078
基因名	TIMP3
Swiss	P35625

别名	TIMP3 (YD35940)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn