

— ABLYBIO, Help Your Research



COL10A1 (YD32489) Rabbit mAb

货号: **AYD16695**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	67kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted, extracellular space, extracellular matrix
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Type X collagen is a product of hypertrophic chondrocytes and has been localized to presumptive mineralization zones of hyaline cartilage
基因ID	3134

基因名	Col10a1
Swiss	Q05306 (https://www.uniprot.org/uniprotkb/Q05306/entry)
别名	COL10A1 (YD32489),COL10A1 (YD32489) Rabbit mAb,Col10a1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)